

WB / Antibody Stripping Buffer (Acid)

Catalog# BWR1076

Size: 250 ml

Lot # Check on the product label

Introduction

WB Stripping Buffer is used for the recycle of deproteinized membrane in WB. In WB, after completing the conjugation of the primary and secondary antibodies as well as the later chemiluminescence, sometimes end user has to detect the proteins with stable expression level like Actin and Tubulin, etc. for reference, or test other proteins for comparison. But for this reagent, through dissociating the conjugation of primary, secondary antibodies and antigens, it removes primary and secondary antibodies from the membrane, then end user can reuse the same membrane to detect other proteins. Thus, compare with running a fresh SDS-PAGE gel, this reagent helps to eliminate the loading error and makes the experiment more comparable and accurate.

For excluding the frequently-used β -Mercaptoethanol, this reagent is nontoxic, tasteless and harmless, can be used at room temperature. This reagent is enough for at least 5 times WB detections on the same membrane, and the whole detection can be completed within 15-30 min.

With the unique component, it is efficient, and has no harm to the proteins.

Kit Components

Components	Size	Storage Instruction
Antibody Stripping Buffer (Acid)	125 ml \times 2	Store at RT for one year.

Protocol

1. Wash membrane with distilled water for 3 times, 5 min each. Discard the distilled water and immerse it into appropriate volume of Antibody Stripping Buffer, incubate at room temperature for 30-60 min and shake slowly.
2. Take out membrane with tweezer. Rinse with TBS-T for 3-4 times, 3-5 min each on shaker.
3. Block with 5% skimmed milk, proceed to the next step.

Notes

1. For horseradish peroxidase system (HRP), 5% skimmed milk should be used as the blocking buffer for each blocking, and for alkaline phosphatase system (AP), casein should be used for each blocking.
2. This Antibody stripping buffer only work on PVDF membrane. The protein transferred to the membrane will has a large loss when use NC membrane.
3. This buffer can work for WB assay by using ECL or the similar chemiluminescence reagents, but not work for non-chemiluminescence reagents, such as DAB, NBT/BCIP.
4. It has slight corrosivity, please wear the lab coat and disposable gloves to operate.

FOR RESEARCH USE ONLY, NOT FOR DIAGNOSTIC AND CLINICAL USE.

Chongqing Biospes Co., Ltd Tel: +86-23-67567091 Fax: +86-23-67745923

7F, Bldg B, High-tech Venture Park, # 107 Erlang Chuangye Rd, Jiulongpo District, Chongqing, 400039, China

www.biospes.com