

Goat anti-human IgG (H+L)-Unconjugated

Catalog# BSA1023	This product was purified from antiserum by immunoaffinity chromatography.
Lot # Check on the product label	
Size 1mg	
Isotype IgG	Specificity It can specifically recognize the heavy chains on human IgG and the light chains on all human immunoglobulins. And Immunoelectrophoresis shows bands only in the gamma region against normal human serum.
Host Goat	
Reactivity Human	Recommend Application Western Blot (1:1,000-5,000) Immunohistochemistry, IHC (1:1,000-5,000) Immunoprecipitation, IP (1:1,000-5,000) Immunofluorescence, IF (1:1,000-5,000) ELISA (1:1,000-5,000) Other applications have not been tested. The optimal dilutions should be determined by end user.
Immunogen Human IgG	
Product Form Liquid	Dilution It is concentrated product, can be diluted with TBST (pH 7.0). Storage Buffer 10mM sodium phosphate, 150mM sodium chloride, pH7.2 Storage Instruction Ship at 4°C, after receipt, aliquot and store at -20°C for long term. Avoid repeated freeze and thaw cycles.

Purification

FOR RESEARCH USE ONLY, NOT FOR DIAGNOSTIC AND CLINICAL USE.